

Product datasheet

Anti-UGT2A3 antibody ab81231

画像数 1

製品の概要

製品名	Anti-UGT2A3 antibody
製品の詳細	Rabbit polyclonal to UGT2A3
アプリケーション	<b>適用あり:</b> WB, ELISA
種交差性	<b>交差種:</b> Human
免疫原	Synthetic peptide designed within residues NVKVILEELI VRGHEVTVLT HSKPSLIDYR KPSALKFEVV HMPQDRTEEN, corresponding to N terminal amino acids 37-86 of Human UGT2A3 (NP_079019) <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	721_B cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab81231** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------


WB	Use a concentration of 1 µg/ml. Predicted molecular weight: 60 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
----	--

アプリケーション	Abreviews	特記事項
----------	-----------	------

ELISA  
 Use at an assay dependent concentration.  
 Titre using peptide based assay: 1/62500.

## ターゲット情報

## 画像



Anti-UGT2A3 antibody (ab81231) at 1 µg/ml  
 (in 5% skim milk / PBS buffer) + 721\_B cell  
 lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000  
 dilution

**Predicted band size** : 60 kDa

**Observed band size** : 60 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors