

Product datasheet

Anti-*Legionella pneumophila* LPS antibody [2F10] ab8263

2 References

製品の概要

製品名	Anti- <i>Legionella pneumophila</i> LPS antibody [2F10]
製品の詳細	Mouse monoclonal [2F10] to <i>Legionella pneumophila</i> LPS
由来種	Mouse
特異性	This antibody reacts with LPS membrane of <i>Legionella pneumophila</i>
アプリケーション	適用あり: ELISA, Sandwich ELISA
種交差性	This antibody does not react with
免疫原	Membrane fraction of <i>Legionella pneumophila</i> , strain Philadelphia 1
特記事項	Detection of <i>Legionella pneumophila</i> .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	PBS with 0.1% sodium azide, pH 7.4
精製度	Protein G purified
一次抗体 備考	Detection of <i>Legionella pneumophila</i> .
ポリモノ	モノクローナル
クローン名	2F10
ミエローマ	Sp2/0
アイソタイプ	IgG1
軽鎖の種類	unknown

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab8263** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		Use at an assay dependent dilution.
Sandwich ELISA		Use at an assay dependent dilution. Can be used as Capture antibody with ab24457 .

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors