

Product datasheet

Rat IgG2a, kappa monoclonal (Spectral Red™) - Isotype Control ab82740

製品の概要

製品名	Rat IgG2a, kappa monoclonal (Spectral Red™) - Isotype Control
製品の詳細	Rat IgG2a, kappa monoclonal (Spectral Red™) - Isotype Control
標識	Spectral Red™
アプリケーション	適用あり: Flow Cyt
特記事項	Conjugated forms should not be frozen and should be protected from prolonged exposure to light.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	Preservative: 0.09% Sodium Azide Constituents: 16% Sucrose, 0.20% Gelatin, PBS and stabilizing agent.
ポリモノ	モノクローナル
クローン名	KLH/G2a-1-1
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab82740** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
Flow Cyt		

追加情報

Flow Cyt: Use 0.2µg for 10⁶ cells.

Not yet tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors