

Product datasheet

Anti-HEATR4 antibody ab83870

画像数 1

製品の概要

製品名	Anti-HEATR4 antibody
製品の詳細	Rabbit polyclonal to HEATR4
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 2 - 51 (VFFSSQYRLH RKSQYLKMAA ANLTFSQEVV WQRGLPSIPY SQYSFDHLYN) of Human HEATR4 (NP_976054).

 [Run BLAST with ExPASy](#)  [Run BLAST with NCBI](#)

ポジティブ・コントロール	293T cell lysate
--------------	------------------

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab83870** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 112 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

関連性	The function of the HEATR4 protein remains unknown.
細胞内局在	Nuclear

画像



Anti-HEATR4 antibody (ab83870) at 1 µg/ml +
293T cell lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 112 kDa

Observed band size: 112 kDa

Additional bands at: < 40 kDa. We are
unsure as to the identity of these extra bands.

Gel concentration 6-18%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors