

## Product datasheet

# Anti-ZNF141 antibody ab85662

### 画像数 1

#### 製品の概要

製品名	Anti-ZNF141 antibody
製品の詳細	Rabbit polyclonal to ZNF141
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide, corresponding to a region within internal sequence amino acids 424-473 (CKECDKAFKQ FSLLSQHKKI HTVDKPYKCK DCDKAFKRFS HLNKHKKIHT) of Human ZNF141 (NP_003432)
ポジティブ・コントロール	human placenta lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: 0.09% Sodium azide Constituents: PBS, 2% Sucrose
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab85662** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 55 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## ターゲット情報

**関連性** A zinc finger encoding cDNA (ZNF141) of the C2-H2/KRAB subfamily has been mapped to the 4p- (Wolf-Hirschhorn) syndrome (WHS) chromosome region. Zinc finger encoding genes would be good candidates for being involved in the multiple developmental defects associated with chromosomal aneusomy—because of their role as transcriptional regulators, their abundance in the genome and their known association with specific developmental disorders.

**細胞内局在** Nuclear

## 画像



Anti-ZNF141 antibody (ab85662) at 1 µg/ml (in 5% skim milk / PBS buffer) + Placenta lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000 dilution

**Predicted band size:** 55 kDa

**Observed band size:** 60 kDa

[why is the actual band size different from the predicted?](#)

**Additional bands at:** 39 kDa, 41 kDa. We are unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors