

## Product datasheet

# Anti-ZNF141 antibody ab85662

### 画像数 1

#### 製品の概要

製品名	Anti-ZNF141 antibody
製品の詳細	Rabbit polyclonal to ZNF141
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide, corresponding to a region within internal sequence amino acids 424-473 (CKECDKAFKQ FSLLSQHKKI HTVDKPYKCK DCDKAFKRFS HLNKHKKIHT) of Human ZNF141 (NP_003432)
ポジティブ・コントロール	human placenta lysate

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab85662** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------


WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 55 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.
----	--	--

## ターゲット情報

**関連性** A zinc finger encoding cDNA (ZNF141) of the C2-H2/KRAB subfamily has been mapped to the 4p- (Wolf-Hirschhorn) syndrome (WHS) chromosome region. Zinc finger encoding genes would be good candidates for being involved in the multiple developmental defects associated with chromosomal aneusomy--because of their role as transcriptional regulators, their abundance in the genome and their known association with specific developmental disorders.

**細胞内局在** Nuclear

## 画像



Anti-ZNF141 antibody (ab85662) at 1 µg/ml  
(in 5% skim milk / PBS buffer) + Placenta  
lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size:** 55 kDa

**Observed band size:** 60 kDa

**Additional bands at:** 39 kDa, 41 kDa. We  
are unsure as to the identity of these extra  
bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors