


Product datasheet

Anti-HOXB2 antibody ab86386

画像数 3

製品の概要

製品名	Anti-HOXB2 antibody
製品の詳細	Rabbit polyclonal to HOXB2
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Horse, Guinea pig, Cow, Cat, Dog 
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 72-121 (AEDGPALPPP PPPPLPAAPP APEFPWMKEK KSAKKPSQSA TSPSPAASAV)of Human HOXB2 (NP_002136). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	HepG2 cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab86386** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
----------	-----------	------

WB Use a concentration of 1 µg/ml. Detects a band of approximately 38 kDa (predicted molecular weight: 38 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

機能	Sequence-specific transcription factor which is part of a developmental regulatory system that provides cells with specific positional identities on the anterior-posterior axis.
配列類似性	Belongs to the Antp homeobox family. Proboscipedia subfamily. Contains 1 homeobox DNA-binding domain.
発生段階	Expressed in whole embryos and fetuses at 5-9 weeks from conception.
細胞内局在	Nucleus.

画像




Anti-HOXB2 antibody (ab86386) at 1 µg/ml +
HepG2 cell lysate at 10 µg


Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 38 kDa



Western Blot using ab86386 at 1.0ug/ml.
Sample type: Human 721_B cell Lysate



Western Blot using ab86386 at 1.0ug/ml.

Sample type: Human Jurkat cell Lysate

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors