

Product datasheet

Anti-FAM81B antibody ab86979

画像数 1

製品の概要

製品名	Anti-FAM81B antibody
製品の詳細	Rabbit polyclonal to FAM81B
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to a region within the N terminal amino acids 2-51 (QLQLGLTLAS SEKRKKSQRL FFKNIKSTKN KAGKASIMSS DTNVNKSASP) of Human FAM81B (NP_689761). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human fetal liver lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab86979** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項


アプリ ケーショ ン	Abreviews	特記事項
------------------	-----------	------

WB Use a concentration of 1 µg/ml. Detects a band of approximately 52 kDa (predicted molecular weight: 52 kDa). Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

配列類似性 Belongs to the FAM81 family.

画像



Anti-FAM81B antibody (ab86979) at 1 µg/ml
(in 5% skim milk / PBS buffer) + Human fetal
liver lysate at 10 µg

Secondary

HRP conjugated anti-Rabbit IgG at 1/50000
dilution

Predicted band size: 52 kDa

Observed band size: 52 kDa

Additional bands at: ~36 kDa. We are
unsure as to the identity of these extra bands.

12% gel

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors