


Product datasheet

Anti-ZNF234 antibody ab87221

画像数 1

製品の概要

|              |                                                                                                                                                                                                                                                             |
|--------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 製品名          | Anti-ZNF234 antibody                                                                                                                                                                                                                                        |
| 製品の詳細        | Rabbit polyclonal to ZNF234                                                                                                                                                                                                                                 |
| アプリケーション     | 適用あり: WB                                                                                                                                                                                                                                                    |
| 種交差性         | 交差種: Human<br>交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Cow, Cat, Dog, Pig                                                                                                 |
| 免疫原          | Synthetic peptide corresponding to a region within internal sequence amino acids 576-625 (GFKWSLNLDM HQRVHTGEKP YTCGECGKHF SQASSLQLHQ SVHTGEKPYK) of Human ZNF234 (NP_006621).<br><a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a> |
| ポジティブ・コントロール | Human fetal pancreas lysate.                                                                                                                                                                                                                                |

製品の特性

|        |                                                                                                |
|--------|------------------------------------------------------------------------------------------------|
| 製品の状態  | Liquid                                                                                         |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. |
| バッファー  | Preservative: None<br>Constituents: 2% Sucrose, PBS                                            |
| 精製度    | Immunogen affinity purified                                                                    |
| ポリモノ   | ポリクローナル                                                                                        |
| アイソタイプ | IgG                                                                                            |

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab87221** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項                                                                                                                                                 |
|----------|-----------|------------------------------------------------------------------------------------------------------------------------------------------------------|
| WB       |           | Use a concentration of 1 µg/ml. Predicted molecular weight: 81 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T. |


## ターゲット情報

---

|       |                                                                                                                                                                                             |
|-------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| 関連性   | ZNF234 belongs to the krueppel C2H2-type zinc-finger protein family. It contains nineteen C2H2-type zinc fingers and one KRAB domain. ZNF234 may be involved in transcriptional regulation. |
| 細胞内局在 | Nuclear                                                                                                                                                                                     |

## 画像

---



Anti-ZNF234 antibody (ab87221) at 1 µg/ml  
(in 5% skim milk / PBS buffer) + Human fetal  
pancreas lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size** : 81 kDa

**Observed band size** : 66 kDa

**Additional bands at** : 40 kDa. We are  
unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors