

## Product datasheet

# Anti-ZNF234 antibody ab87221

画像数 1

### 製品の概要

製品名	Anti-ZNF234 antibody
製品の詳細	Rabbit polyclonal to ZNF234
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Cow, Cat, Dog, Pig
免疫原	Synthetic peptide corresponding to a region within internal sequence amino acids 576-625 (GFKWLSNLDLDM HQRVHTGEKP YTCGECGKHF SQASSLQLHQ SVHTGEKPYK) of Human ZNF234 (NP_006621). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	Human fetal pancreas lysate.

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab87221** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 81 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.


## ターゲット情報

---

関連性	ZNF234 belongs to the krueppel C2H2-type zinc-finger protein family. It contains nineteen C2H2-type zinc fingers and one KRAB domain. ZNF234 may be involved in transcriptional regulation.
細胞内局在	Nuclear

## 画像

---



Anti-ZNF234 antibody (ab87221) at 1 µg/ml  
(in 5% skim milk / PBS buffer) + Human fetal  
pancreas lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size** : 81 kDa

**Observed band size** : 66 kDa

**Additional bands at** : 40 kDa. We are  
unsure as to the identity of these extra bands.

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors