

Product datasheet

Anti-EN1 / Engrailed 1 antibody ab87963

画像数 2

製品の概要

製品名	Anti-EN1 / Engrailed 1 antibody
製品の詳細	Mouse monoclonal to EN1 / Engrailed 1
由来種	Mouse
アプリケーション	適用あり: WB, ELISA
種交差性	交差種: Human
免疫原	Recombinant fragment, corresponding to amino acids 266-392 of Human EN1/ Engrailed 1 (NP_001417) with a 26 kDa tag
ポジティブ・コントロール	Jurkat cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	モノクローナル
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab87963** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 40 kDa.
ELISA		Use at an assay dependent dilution.

ターゲット情報


配列類似性

Belongs to the engrailed homeobox family.
Contains 1 homeobox DNA-binding domain.

細胞内局在

Nucleus.

画像




Anti-EN1 / Engrailed 1 antibody (ab87963) at
1 µg/ml + Immunogen (Recombinant protein)
at 0.1 µg

Developed using the ECL technique.

Predicted band size: 40 kDa

Observed band size: 40 kDa



Anti-EN1 / Engrailed 1 antibody (ab87963) at
1 µg/ml + Jurkat cell lysate at 50 µg

Developed using the ECL technique.

Predicted band size: 40 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Promise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors