

Product datasheet

Anti-RNF23 antibody ab88707

画像数 3

製品の概要

製品名	Anti-RNF23 antibody
製品の詳細	Mouse polyclonal to RNF23
アプリケーション	適用あり: WB, ICC/IF
種交差性	交差種: Human
免疫原	Full length Human RNF23 (AAH07661.1)
ポジティブ・コントロール	Human kidney tissue lysate, RNF23 transfected 293T cell lysate and HeLa cells.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab88707** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Detects a band of approximately 54 kDa (predicted molecular weight: 54 kDa).
ICC/IF		Use a concentration of 10 µg/ml.

ターゲット情報


機能	May facilitate apoptosis by inhibiting the poly-ubiquitination and subsequent proteasome-mediated degradation of the pro-apoptotic protein MOAP1.
組織特異性	Ubiquitous; highly expressed in brain, heart, kidney, liver, skeletal muscle, spleen and testis.
配列類似性	<p>Belongs to the TRIM/RBCC family.</p> <p>Contains 1 B box-type zinc finger.</p> <p>Contains 1 B30.2/SPRY domain.</p> <p>Contains 1 RING-type zinc finger.</p>
細胞内局在	Cytoplasm > cytosol. Mitochondrion. Found predominantly in the cytosol. Partial shift from the cytosol to the mitochondria when co-localized with MOAP1.

Anti-RNF23 antibody 画像



Anti-RNF23 antibody (ab88707) at 1 µg/ml +
Human kidney tissue lysate at 50 µg

Predicted band size : 54 kDa
Observed band size : 54 kDa




All lanes : Anti-RNF23 antibody (ab88707) at
1 µg/ml

Lane 1 : RNF23 transfected 293T cell lysate
Lane 2 : Non-transfected lysate

Lysates/proteins at 25 µg per lane.

Predicted band size : 54 kDa
Observed band size : 54 kDa
Additional bands at : 33 kDa, 35 kDa. We
are unsure as to the identity of these extra
bands.



ab88707 at 10μg/ml staining RNF23 in HeLa cells by Immunofluorescence.

Immunocytochemistry/ Immunofluorescence -
RNF23 antibody (ab88707)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors