

Product datasheet

Recombinant Human RPLP2 protein ab88711

製品の概要

製品名	Recombinant Human RPLP2 protein
タンパク質長	Full length protein

製品の詳細

由来	Recombinant
由来	Saccharomyces cerevisiae

アミノ酸配列

生物種	Human
配列	MRYVASYLLAALGGNSSPSAKDIKKILDSVGIEADDDRLNKVISELNGKN IEDVIAQGIGKLASVPAGGAVAVSAAPGSAAPAAGSAPAAAEEKKDEKKE ESEESDDDMGFGLFD

特性

Our [Abpromise guarantee](#) covers the use of **ab88711** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	SDS-PAGE
精製度	> 90 % SDS-PAGE. ab88711 is 90% pure as determined by SDS-PAGE.
製品の状態	Liquid

前処理および保存

保存方法および安定性	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles. Preservative: None Constituents: 30% Glycerol, 0.5% Triton-X-100, 50mM HEPES, 30mM Glutathione, 100mM Sodium chloride, 1mM DTT, pH 7.5
------------	---

関連情報

機能	Plays an important role in the elongation step of protein synthesis.
配列類似性	Belongs to the ribosomal protein L12P family.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors