

Product datasheet

Anti-RBM46 antibody ab88733

画像数 2

製品の概要

製品名	Anti-RBM46 antibody
製品の詳細	Mouse polyclonal to RBM46
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length Human RBM46 protein (NP_659416).
ポジティブ・コントロール	Human pancreas tissue lysate; RBM46 transfected 293T cell lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab88733** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Detects a band of approximately 60 kDa (predicted molecular weight: 60 kDa).

ターゲット情報

関連性

Contains 3 RRM (RNA recognition motif) domains. Molecular function: RNA binding.


画像



Anti-RBM46 antibody (ab88733) at 1/500 dilution + Human pancreas tissue lysate at 50 μ g

Predicted band size: 60 kDa

Observed band size: 60 kDa



All lanes : Anti-RBM46 antibody (ab88733) at 1/500 dilution

Lane 1 : RBM46 transfected 293T cell lysate

Lane 2 : Non transfected 293T cell lysate

Lysates/proteins at 25 μ g per lane.

Predicted band size: 60 kDa

Observed band size: 60 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors