

Product datasheet

Anti-Surf6 antibody ab88832

画像数 2

製品の概要

製品名	Anti-Surf6 antibody
製品の詳細	Mouse polyclonal to Surf6
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length Human Surf6 (NP_006744, amino acids 1-361)
ポジティブ・コントロール	Human liver tissue lysate, Surf6 transfected 293T cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: None Constituents: Whole serum
精製度	Whole antiserum
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab88832** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 41 kDa.

ターゲット情報

機能	Binds to both DNA and RNA in vitro, with a stronger binding capacity for RNA. May represent a nucleolar constitutive protein involved in ribosomal biosynthesis or assembly.
----	--


## 配列類似性

Belongs to the SURF6 family.

## 細胞内局在

Nucleus > nucleoplasm. Nucleus > nucleolus. Granular component of the nucleolus.

## 画像




Western blot - Surf6 antibody (ab88832)

Anti-Surf6 antibody (ab88832) at 1/500 dilution + human liver tissue lysate at 50 µg

**Predicted band size** : 41 kDa

**Observed band size** : 41 kDa



Western blot - Surf6 antibody (ab88832)

**All lanes** : Anti-Surf6 antibody (ab88832) at 1/500 dilution

**Lane 1** : Surf6 transfected 293T cell lysate

**Lane 2** : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

**Predicted band size** : 41 kDa

**Observed band size** : 45 kDa

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors