

Product datasheet

Anti-HLA Class II DRB1 antibody ab89315

1 References [画像数 2](#)

製品の概要

製品名	Anti-HLA Class II DRB1 antibody
製品の詳細	Mouse polyclonal to HLA Class II DRB1
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Full length protein corresponding to Human HLA Class II DRB1 aa 1-266. Database link: P01911
ポジティブ・コントロール	Human spleen tissue lysate and HLA DRB1 transfected 293T cell lysate.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 1X PBS, pH 7.2
精製度	Protein A purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


Our [Abpromise guarantee](#) covers the use of **ab89315** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 30 kDa.

ターゲット情報

画像




Western blot - Anti-HLA Class II DRB1 antibody (ab89315)

Anti-HLA Class II DRB1 antibody (ab89315)
at 1 µg/ml + Human spleen tissue lysate at 50 µg

Predicted band size: 30 kDa

Observed band size: 30 kDa



Western blot - Anti-HLA Class II DRB1 antibody (ab89315)

All lanes : Anti-HLA Class II DRB1 antibody (ab89315) at 1 µg/ml

Lane 1 : HLA DRB1 transfected 293T cell lysate

Lane 2 : Non-transfected 293T cell lysate

Lysates/proteins at 25 µg per lane.

Predicted band size: 30 kDa

Observed band size: 30 kDa

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors