

Product datasheet

Anti-Macrophage Specific Gene antibody ab90019

画像数 1

製品の概要

製品名	Anti-Macrophage Specific Gene antibody
製品の詳細	Rabbit polyclonal to Macrophage Specific Gene
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Synthetic peptide corresponding to a region within amino acids 667-716 (ILAVVITLAI YGTRKFKKKA YQAI EERQSL VPGTAATGDT TYQEQQSPA) of Human Macrophage Specific Gene (NP_001034485). Run BLAST with ExPASy Run BLAST with NCBI
ポジティブ・コントロール	Human fetal spleen lysate

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
特記事項 (精製)	ab90019 is purified by a peptide affinity chromatography method.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab90019** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.


アプリケーション	Abreviews	特記事項
----------	-----------	------

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 78 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

ターゲット情報

組織特異性	Expressed in macrophages and peripheral blood monocytes.
配列類似性	Contains 1 MACPF domain.
細胞内局在	Membrane.

画像



Western blot - Macrophage Specific Gene antibody (ab90019)

Anti-Macrophage Specific Gene antibody (ab90019) at 1 µg/ml (in 5% skim milk / PBS buffer) + Human fetal spleen lysate at 10 µg

Secondary

anti-Rabbit IgG HRP at 1/50000 dilution

Predicted band size: 78 kDa

Gel concentration: 12%

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors