

## Product datasheet

## Anti-SEMA6C antibody [MM0536-9S21] ab90240

## 製品の概要

製品名	Anti-SEMA6C antibody [MM0536-9S21]
製品の詳細	Mouse monoclonal [MM0536-9S21] to SEMA6C
由来種	Mouse
アプリケーション	適用あり: WB
種交差性	交差種: Human
免疫原	Recombinant fragment corresponding to the extracellular domain of Human SEMA6C

## 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	Constituent: PBS
精製度	Protein G purified
特記事項(精製)	The IgG fraction of culture supernatant was purified by Protein G affinity chromatography and filtered through a 0.2 µm filter.
ポリ/モノ	モノクローナル
クローン名	MM0536-9S21
アイソタイプ	IgG2

## アプリケーション

## The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab90240の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/1000. Predicted molecular weight: 100 kDa.

## ターゲット情報

## 関連性

SEMA6C belongs to a subfamily characterized by an extracellular semaphorin domain, a transmembrane domain, and a long cytoplasmic tail. Members of this class can repel sympathetic and dorsal root ganglion axons in vitro, consistent with a traditional role as guidance signals. SEMA6C shows growth cone collapsing activity on dorsal root ganglion (DRG) neurons in vitro and may be a stop signal for the DRG neurons in their target areas, and possibly also for other neurons. SEMA6C may also be involved in the maintenance and remodeling of neuronal connections.

## 細胞内局在

Membrane; Single-pass type I membrane protein.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

---

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

---

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors