


## Product datasheet

# Anti-PABPC5 antibody ab94568

画像数 1

### 製品の概要

製品名	Anti-PABPC5 antibody
製品の詳細	Rabbit polyclonal to PABPC5
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit, Horse, Guinea pig, Cow, Cat, Dog, Chimpanzee 
免疫原	Synthetic peptide corresponding to a region within N terminal amino acids 72-121 (AEWALNTMNF DLINGKPFRL MWSQPDDRLR KSGVGNIFIK NLDKSIDNRA) of Human PABPC5 (NP_543022). <a href="#">Run BLAST with ExPASy</a> <a href="#">Run BLAST with NCBI</a>
ポジティブ・コントロール	NCH226 cell lysate

### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid repeated freeze / thaw cycles.
バッファー	Preservative: None Constituents: 2% Sucrose, PBS
精製度	Immunogen affinity purified
ポリモノ	ポリクローナル
アイソタイプ	IgG

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab94568** in the following tested applications.


The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		Use a concentration of 1 µg/ml. Predicted molecular weight: 43 kDa. Good results were obtained when blocked with 5% non-fat dry milk in 0.05% PBS-T.

## ターゲット情報

関連性	PABPC5 (polyadenylate-binding protein 5) binds the polyadenylate tail of mRNA and is thought to be involved in cytoplasmic regulatory processes of mRNA metabolism. In vivo, it may also bind to cytoplasmic RNA sequences other than polyadenylate.
細胞内局在	Cytoplasmic

## 画像



Anti-PABPC5 antibody (ab94568) at 1 µg/ml  
(in 5% skim milk / PBS buffer) + NCI-H226 cell  
lysate at 10 µg

### Secondary

HRP conjugated anti-Rabbit IgG at 1/50000  
dilution

**Predicted band size: 43 kDa**

Gel concentration: 12%

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors