

Product datasheet

Anti-HLA G antibody [MEM-G/9] (FITC) ab94890

1 References

製品の概要

製品名	Anti-HLA G antibody [MEM-G/9] (FITC)
製品の詳細	Mouse monoclonal [MEM-G/9] to HLA G (FITC)
由来種	Mouse
標識	FITC. Ex: 493nm, Em: 528nm
特異性	The antibody reacts with the native form of Human HLA G1 on the cell surface as well as with soluble HLA G5 isoform in its beta2-microglobulin associated form.
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Human
免疫原	Recombinant Human HLA G refolded with beta2-microglobulin and peptide.
特記事項	The purified antibody is conjugated with Fluorescein isothiocyanate (FITC) under optimum conditions. The reagent is free of unconjugated FITC.

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 15mM Sodium Azide Constituents: PBS, pH 7.4
精製度	Ascites
特記事項(精製)	Protein G purified
ポリ/モノ	モノクローナル
クローン名	MEM-G/9
アイソタイプ	IgG1

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab94890** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
----------	-----------	------

Flow Cyt

1/500.

[ab91356](#) - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.

ターゲット情報

機能	Involved in the presentation of foreign antigens to the immune system. Plays a role in maternal tolerance of the fetus by mediating protection from the deleterious effects of natural killer cells, cytotoxic T lymphocytes, macrophages and mononuclear cells.
組織特異性	Expressed in trophoblasts.
配列類似性	Belongs to the MHC class I family. Contains 1 Ig-like C1-type (immunoglobulin-like) domain.
細胞内局在	Membrane.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors