

## Product datasheet

# Anti-TRIM antibody [TRIM-04] ab9529

### 3 References

#### 製品の概要

製品名	Anti-TRIM antibody [TRIM-04]
製品の詳細	Mouse monoclonal [TRIM-04] to TRIM
由来種	Mouse
特異性	Human T cell receptor-interacting molecule (TRIM)
アプリケーション	<b>適用あり:</b> WB, IP
種交差性	<b>交差種:</b> Human
免疫原	Recombinant fragment, corresponding to amino acids 29-186 of Human TRIM.

#### 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 15mM Sodium Azide Constituents: PBS, pH 7.4
精製度	>95% by SDS-PAGE
特記事項 (精製)	Purified from ascites by protein-A affinity chromatography.
ポリ/モノ	モノクローナル
クローン名	TRIM-04
アイソタイプ	IgG2a

#### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab9529** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
WB		

IP

**追加情報**

IP: Use at an assay dependant concentration.

WB: Use at an assay dependant concentration. Use under non-reducing conditions.

Not tested in other applications.

Optimal dilutions/concentrations should be determined by the end user.

**ターゲット情報****機能** Stabilizes the TCR (T-cell antigen receptor)/CD3 complex at the surface of T-cells.**組織特異性** Strongly expressed in thymus, and to a lesser extent in spleen, lymph node and peripheral blood lymphocytes. Present in T-cells and NK cells, but not B-cells (at protein level).**発生段階** Strongly expressed in fetal thymus at weeks 17-24 of gestation. Undetectable in bone marrow and fetal liver.**翻訳後修飾** Phosphorylated on tyrosines by LCK or FYN upon TCR activation.**細胞内局在** Cell membrane.**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"**Our Abpromise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors