

Product datasheet

Mouse monoclonal [G1 1C5] Anti-Rat IgG1 heavy chain (FITC) ab99659

製品の概要

製品名	Mouse monoclonal [G1 1C5] Anti-Rat IgG1 heavy chain (FITC)
由来種	Mouse
ターゲット生物種	Rat
特異性	Reacts with the gamma 1 heavy chain (Fc) of Rat IgG1.
アプリケーション	適用あり: ICC/IF, Flow Cyt, ELISA
免疫原	The details of the immunogen for this antibody are not available.
標識	FITC. Ex: 493nm, Em: 528nm

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: 0.1% Sodium Azide Constituents: PBS
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	G1 1C5
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99659** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 0.3µg for 10 ⁶ cells. ab106163 - Mouse monoclonal IgG1, is suitable for use as an isotype control with this antibody.
ELISA		Use at an assay dependent concentration.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors