

Product datasheet

Mouse monoclonal K4F5 Anti-Rat kappa light chain (HRP)
ab99692

画像数 1

製品の概要

製品名	Mouse monoclonal K4F5 Anti-Rat kappa light chain (HRP)
由来種	Mouse
ターゲット生物種	Rat
アプリケーション	適用あり: ELISA, WB
免疫原	The details of the immunogen for this antibody are not available.
標識	HRP

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C.
バッファー	Preservative: None Constituents: 50% PBS, 50% Glycerol, pH 7.4
精製度	Protein G purified
ポリモノ	モノクローナル
クローン名	K4F5
アイソタイプ	IgG1
軽鎖の種類	kappa


アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99692** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ELISA		1/2000 - 1/5000.
WB		Use at an assay dependent concentration.

画像



Lane 1: Rat IgG2a Kappa loaded at 0.2ug/well

Lane 2: Rat IgM Lambda loaded at 0.2ug/well
ab99692 diluted to 1:2,000

Western blot - Mouse monoclonal [K4F5] Secondary
Antibody to Rat kappa - light chain (HRP) (ab99692)

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors