

## Product datasheet

# Mouse monoclonal A-1 Anti-Chicken IgA heavy chain ab99726

### 製品の概要

|          |                                                           |
|----------|-----------------------------------------------------------|
| 製品名      | Mouse monoclonal A-1 Anti-Chicken IgA heavy chain         |
| 由来種      | Mouse                                                     |
| ターゲット生物種 | Chicken                                                   |
| 特異性      | ab99726 reacts with the alpha heavy chain of Chicken IgA. |
| アプリケーション | <b>適用あり:</b> ICC/IF, IHC-Fr, ELISA, Flow Cyt              |

### 製品の特性

|        |                                                                                                                                  |
|--------|----------------------------------------------------------------------------------------------------------------------------------|
| 製品の状態  | Liquid                                                                                                                           |
| 保存方法   | Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle. |
| バッファー  | pH: 8.20<br>Constituent: 100% Borate buffered saline<br><br>No preservatives or amine-containing buffer salts added.             |
| 精製度    | IgG fraction                                                                                                                     |
| ポリ/モノ  | モノクローナル                                                                                                                          |
| クローン名  | A-1                                                                                                                              |
| アイソタイプ | IgG2b                                                                                                                            |
| 軽鎖の種類  | kappa                                                                                                                            |

### アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99726** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

| アプリケーション | Abreviews | 特記事項                                                            |
|----------|-----------|-----------------------------------------------------------------|
| ICC/IF   |           | Use at an assay dependent concentration.                        |
| IHC-Fr   |           | Use at an assay dependent concentration. Acetone fixed sections |

| アプリケーション | Abreviews | 特記事項                                                                                                                                                     |
|----------|-----------|----------------------------------------------------------------------------------------------------------------------------------------------------------|
| ELISA    |           | Use at an assay dependent concentration.                                                                                                                 |
| Flow Cyt |           | Use at an assay dependent concentration. <a href="#">ab170192</a> -Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody. |

**Please note:** All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

### Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

### Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors