

Product datasheet

Mouse monoclonal [A-1] Anti-Chicken IgA heavy chain ab99726

製品の概要

製品名	Mouse monoclonal [A-1] Anti-Chicken IgA heavy chain
由来種	Mouse
ターゲット生物種	Chicken
特異性	ab99726 reacts with the alpha heavy chain of Chicken IgA.
アプリケーション	適用あり: ICC/IF, IHC-Fr, ELISA, Flow Cyt

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	pH: 8.20 Constituent: 100% Borate buffered saline No preservatives or amine-containing buffer salts added.
精製度	IgG fraction
ポリ/モノ	モノクローナル
クローン名	A-1
アイソタイプ	IgG2b
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99726** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC/IF		Use at an assay dependent concentration.
IHC-Fr		Use at an assay dependent concentration. Acetone fixed sections
ELISA		Use at an assay dependent concentration.

Flow Cyt

Use at an assay dependent concentration.

[ab170192](#) - Mouse monoclonal IgG2b, is suitable for use as an isotype control with this antibody.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors