

Product datasheet

Mouse monoclonal HP6029 Anti-Human IgE Fc (Alkaline Phosphatase) ab99835

製品の概要

製品名	Mouse monoclonal HP6029 Anti-Human IgE Fc (Alkaline Phosphatase)
ターゲット生物種	Human
特異性	Reacts with the Fc portion of the heavy chain of human IgE as demonstrated by ELISA
アプリケーション	適用あり: ICC, Dot blot, ELISA, WB, IHC-Fr
免疫原	Other Immunogen Type corresponding to IgE.
標識	Alkaline Phosphatase

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.1% Sodium Azide Constituents: 50% Glycerol, 50mM Tris, 1mM Magnesium chloride, pH 8.0
精製度	Protein A purified
特記事項 (精製)	Source: Ascites fluid
ポリ/モノ	モノクローナル
クローン名	HP6029
アイソタイプ	IgG2a
軽鎖の種類	kappa

アプリケーション

Our [Abpromise guarantee](#) covers the use of **ab99835** in the following tested applications.

The application notes include recommended starting dilutions; optimal dilutions/concentrations should be determined by the end user.

アプリケーション	Abreviews	特記事項
ICC		Use at an assay dependent dilution.

アプリケーション	Abreviews	特記事項
Dot blot		Use at an assay dependent dilution.
ELISA		1/500 - 1/1000.
WB		Use at an assay dependent dilution.
IHC-Fr		Use at an assay dependent dilution.

Please note: All products are "FOR RESEARCH USE ONLY AND ARE NOT INTENDED FOR DIAGNOSTIC OR THERAPEUTIC USE"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <http://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors