

Product datasheet

LY 341495 disodium salt (mM/ml), group II/group III mGlu antagonist (water-soluble) ab144500

画像数 1

製品の概要


製品名	LY 341495 disodium salt (mM/ml), group II/group III mGlu antagonist (water-soluble)
製品の詳細	Water-soluble group II/group III mGlu antagonist. Disodium salt of LY 341495 (ab120199). 1 ml water soluble pack.
生理活性の詳細	Very selective nanomolar potent antagonist for group II mGlu receptors. Also a relatively potent antagonist for group III mGlu receptors at high nanomolar to low micromolar concentrations. Disodium salt of LY 341495 (ab120199) which displays K_d values of 1.67, 0.75, 31.6, 72.7, 14 nM at human mGlu ₂ , mGlu ₃ , mGlu ₆ , mGlu _{7a} , and human mGlu _{8a} receptors respectively. LY 341495 displays IC ₅₀ values of 7.8, 8.2 and 22 μ M for mGlu _{1a} , mGlu _{5a} and mGlu ₄ receptors. It is readily brain penetrant and active <i>in vivo</i> .

Soluble in 1 ml water to give specified mM/ml concentration. [Find out more.](#)

精製度

> 98%

構造式



製品の特性

体系名	(2S)-2-Amino-2-[(1S,2S)-2-carboxycycloprop-1-yl]-3-(xanth-9-yl)propanoic acid disodium salt
分子量	397.34
分子式	C ₂₀ H ₁₇ NNa ₂ O ₅
保存方法	Store at -20°C. Store under desiccating conditions. The product can be stored for up to 12 months.
溶解性	Soluble in 1 ml water to give specified mM/ml concentration
使用に関する注意	Wherever possible, you should prepare and use solutions on the same day. However, if you need to make up stock solutions in advance, we recommend that you store the solution as aliquots in tightly sealed vials at -20°C. Generally, these will be useable for up to one month. Before use, and prior to opening the vial we recommend that you allow your product to equilibrate to room


temperature for at least 1 hour.

Need more advice on solubility, usage and handling? Please visit our [frequently asked questions \(FAQ\) page](#) for more details.

由来

Synthetic

画像



2D chemical structure image of ab144500, LY 341495 disodium salt (mM/ml), group II/group III mGlu antagonist (water-soluble)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES, NOT FOR USE IN HUMANS"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors
- Abcam biochemicals are novel compounds and we have not tested their biological activity in house. Please use the literature to identify how to use these products effectively. If you require further assistance please contact the scientific support team