

Streptavidin-HRP Solution ab210901

3 References

製品の概要

製品名	Streptavidin-HRP Solution
製品の概要	<p>Streptavidin-HRP Solution is recommended for use with the Matched Antibody Pair kits for detecting biotinylated molecules during ELISA development.</p> <p>See our website for a protocol.</p> <p>We recommend an initial 1:500 - 1:5,000 dilution range of the Streptavidin-HRP Solution. Optimal dilution should however be determined by the end-user.</p> <p>The Streptavidin-HRP Solution is also available as part of our Matched Antibody Pair Kit Accessory Pack (ab210905).</p> <p>Concentration: 50 µg/mL</p>

製品の特性

保存方法	Store at -20°C. Please refer to protocols.
------	--

内容	110 µl
Streptavidin HRP Solution (ab210901) 50 µg/mL	1 x 110µl

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet

- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors