

Bcl G/BCL2L14 overexpression 293T lysate (whole cell) ab94328

画像数 2

医薬用外劇物

製品の概要

製品名	Bcl G/BCL2L14 overexpression 293T lysate (whole cell)
特記事項	ab94328 is a 293T cell transfected lysate in which Human Bcl G/BCL2L14 has been transiently over-expressed using a pCMV-Bcl G/BCL2L14 plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB

製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Bcl-G is a novel pro-apoptotic member of the Bcl-2 family. The human Bcl-G gene encodes two proteins through alternative mRNA splicing, Bcl-G(L) (long) and Bcl-G(S) (short) consisting of 327 and 252 amino acids in length, respectively. Bcl-G(L) mRNA is expressed widely in adult human tissues, whereas Bcl-G(S) mRNA was found only in testis. Overexpression of Bcl-G(L) or Bcl-G(S) in cells induced apoptosis although Bcl-G(S) was far more potent than Bcl-G(L).

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab94328の使用に適用されます
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

画像




SDS-PAGE -

Bcl G/BCL2L14 overexpression 293T lysate (whole cell)

(ab94328)

ab94328 at 15µg/lane on an SDS-PAGE gel.



Western blot -

Bcl G/BCL2L14 overexpression 293T lysate (whole cell)

(ab94328)

All lanes : Anti-Bcl G/BCL2L14 antibody (**ab57742**) at 1/500 dilution

Lane 1 : Bcl G/BCL2L14 overexpression 293T lysate (whole cell) (ab94328)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors