

### Cathepsin H overexpression 293T lysate (whole cell) ab94087

画像数 2

医薬用外劇物

#### 製品の概要

製品名	Cathepsin H overexpression 293T lysate (whole cell)
特記事項	ab94087 is a 293T cell transfected lysate in which Human Cathepsin H has been transiently overexpressed using a pCMV-Cathepsin H plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB

#### 製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Important for the overall degradation of proteins in lysosomes. Similarity: Belongs to the peptidase C1 family.

#### アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab94087の使用に適用されず**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

#### 画像



SDS-PAGE - Cathepsin H overexpression 293T lysate (whole cell) (ab94087)

ab94087 at 15µg/lane on an SDS-PAGE gel.



Western blot - Cathepsin H overexpression 293T lysate (whole cell) (ab94087)

**All lanes** : Anti-Cathepsin H antibody ([ab88683](#)) at 1/500 dilution

**Lane 1** : Cathepsin H overexpression 293T lysate (whole cell) (ab94087)

**Lane 2** : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

**Secondary**

**All lanes** : Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

**Our Promise to you: Quality guaranteed and expert technical support**

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

**Terms and conditions**

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors