

Cytohesin 2 overexpression 293T lysate (whole cell) ab94253

画像数 2

医薬用外劇物

製品の概要

製品名	Cytohesin 2 overexpression 293T lysate (whole cell)
特記事項	ab94253 is a 293T cell transfected lysate in which Human Cytohesin 2 has been transiently over-expressed using a pCMV-Cytohesin 2 plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB

製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: Acts as a guanine-nucleotide exchange factor (GEF). Promotes guanine-nucleotide exchange on ARF1, ARF3 and ARF6. Promotes the activation of ARF factors through replacement of GDP with GTP. The cell membrane form, in association with ARL4 proteins, recruits ARF6 to the plasma membrane. Tissue specificity: Ubiquitous. Similarity: Contains 1 PH domain. Contains 1 SEC7 domain. Domain: The PH domain is necessary and sufficient for plasma membrane relocalization.

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab94253の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

画像



SDS-PAGE - Cytohesin 2 overexpression 293T
lysate (whole cell) (ab94253)

ab94253 at 15µg/lane on an SDS-PAGE gel



Western blot - Cytohesin 2 overexpression 293T
lysate (whole cell) (ab94253)

All lanes : Anti-Cytohesin 2 antibody (**ab56510**) at 1/500 dilution

Lane 1 : Cytohesin 2 overexpression 293T lysate (whole cell)
(ab94253)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at
1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors