

### Flotillin 2 / ESA overexpression 293T lysate (whole cell) ab94260

画像数 2

医薬用外劇物

#### 製品の概要

製品名	Flotillin 2 / ESA overexpression 293T lysate (whole cell)
特記事項	ab94260 is a 293T cell transfected lysate in which Human Flotillin 2 / ESA has been transiently overexpressed using a pCMV-Flotillin 2 / ESA plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB

#### 製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: May act as a scaffolding protein within caveolar membranes, functionally participating in formation of caveolae or caveolae-like vesicles. May be involved in epidermal cell adhesion and epidermal structure and function. Tissue specificity: In skin, expressed in epidermis and epidermal appendages but not in dermis. Expressed in all layers of the epidermis except the basal layer. In hair follicles, expressed in the suprabasal layer but not the basal layer. Also expressed in melanoma and carcinoma cell lines, fibroblasts and foreskin melanocytes. Similarity: Belongs to the band 7/mec-2 family. Flotillin subfamily.

#### アプリケーション

**The Abpromise guarantee**      **Abpromise保証は、次のテスト済みアプリケーションにおけるab94260の使用に適用されず**  
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

## 画像



SDS-PAGE - Flotillin 2 / ESA overexpression 293T lysate (whole cell) (ab94260)

ab94260 at 15µg/lane on an SDS-PAGE gel.



Western blot - Flotillin 2 / ESA overexpression 293T lysate (whole cell) (ab94260)

**All lanes :** Anti-Flotillin 2/ESA antibody (**ab88498**) at 1/500 dilution

**Lane 1 :** Flotillin 2 / ESA overexpression 293T lysate (whole cell) (ab94260)

**Lane 2 :** 293T Non Transfected Lysate

Lysates/proteins at 25 µg per lane.

### Secondary

**All lanes :** Goat Anti-mouse IgG (H and L) HRP conjugated at 1/2500 dilution

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors