

Granzyme B overexpression 293T lysate (whole cell) ab94040

画像数 2

医薬用外劇物

製品の概要

製品名	Granzyme B overexpression 293T lysate (whole cell)
特記事項	ab94040 is a 293T cell transfected lysate in which Human Granzyme B has been transiently over-expressed using a pCMV-Granzyme B plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB

製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Function: This enzyme is necessary for target cell lysis in cell-mediated immune responses. It cleaves after Asp. Seems to be linked to an activation cascade of caspases (aspartate-specific cysteine proteases) responsible for apoptosis execution. Cleaves caspase-3, -7, -9 and 10 to give rise to active enzymes mediating apoptosis. Similarity: Belongs to the peptidase S1 family. Granzyme subfamily. Contains 1 peptidase S1 domain.

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab94040の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

画像



SDS-PAGE - Granzyme B overexpression 293T
lysate (whole cell) (ab94040)

ab94040 at 15µg/lane on an SDS-PAGE gel.



Western blot - Granzyme B overexpression 293T
lysate (whole cell) (ab94040)

All lanes : Anti-Granzyme B antibody (**ab89415**) at 1/500 dilution

Lane 1 : Granzyme B 293T Transfected Lysate - Positive Control
(ab94040)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Anti-mouse IgG (H and L) HRP conjugated at 1/2500
dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors