

Liver Arginase overexpression 293T lysate (whole cell) ab94083

画像数 2

医薬用外劇物

製品の概要

製品名	Liver Arginase overexpression 293T lysate (whole cell)
特記事項	ab94083 is a 293T cell transfected lysate in which Human liver Arginase has been transiently over-expressed using a pCMV-ARG1 Full length plasmid. The lysate is provided in 1 x Sample Buffer. Note: For more detailed how the transfected lysate was prepared view preparation notes .
アプリケーション	適用あり: WB

製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Pathway: Nitrogen metabolism; urea cycle; L-ornithine and urea from L-arginine: step 1/1. Disease: Defects in ARG1 are the cause of argininemia (ARGIN) [MIM:207800]; also known as hyperargininemia. Argininemia is a rare autosomal recessive disorder of the urea cycle. Arginine is elevated in the blood and cerebrospinal fluid, and periodic hyperammonemia occurs. Clinical manifestations include developmental delay, seizures, mental retardation, hypotonia, ataxia, progressive spastic quadriplegia. Similarity: Belongs to the arginase family.

アプリケーション

The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab94083の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

画像



SDS-PAGE - Liver Arginase overexpression 293T
lysate (whole cell) (ab94083)

ab94083 at 15µg/lane on an SDS-PAGE gel.



Western blot - Liver Arginase overexpression 293T
lysate (whole cell) (ab94083)

All lanes : Anti-liver Arginase antibody ([ab89184](#)) at 1/500 dilution

Lane 1 : Liver Arginase overexpression 293T lysate (whole cell)
(ab94083)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at
1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors