

p19 INK4d overexpression 293T lysate (whole cell) ab94315

画像数 2

医薬用外劇物

製品の概要

製品名	p19 INK4d overexpression 293T lysate (whole cell)
特記事項	ab94315 is a 293T cell transfected lysate in which Human p19 INK4d has been transiently over-expressed using a pCMV-p19 INK4d plasmid. The lysate is provided in 1X Sample Buffer.
アプリケーション	適用あり: WB

製品の特性


Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	Constituents: 0.01% Bromophenol blue, 2.3% Beta mercaptoethanol, 2% Sodium lauryl sulfate, 0.788% Tris HCl, 10% Glycerol (glycerin, glycerine)
背景	Cyclins are important in regulating the cell cycle through their formation of enzymatic complexes with various cyclin dependent kinases. The activity of the cyclin D dependent kinases is controlled by inhibitors such as the INK4 family, which includes INK4a, INK4b, INK4c, and INK4d (p19 INK4d). p19 INK4d is a specific inhibitor of cdk 4 / cdk 6 and is periodically expressed during cell cycle. It is maximally induced as cells enter S phase.

アプリケーション

The Abpromise guarantee Abpromise保証は、 次のテスト済みアプリケーションにおけるab94315の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。


アプリケーション	Abreviews	特記事項
WB		Use at an assay dependent dilution.

画像



SDS-PAGE - p19 INK4d overexpression 293T lysate
(whole cell) (ab94315)

ab94315 at 15µg/lane on an SDS-PAGE gel



Western blot - p19 INK4d overexpression 293T
lysate (whole cell) (ab94315)

All lanes : Anti-p19 INK4d antibody (**ab77223**) at 1/500 dilution

Lane 1 : p19 INK4d overexpression 293T lysate (whole cell)
(ab94315)

Lane 2 : 293T non-transfected lysate

Lysates/proteins at 25 µg per lane.

Secondary

All lanes : Goat Anti-mouse IgG (H and L) HRP conjugated at
1/2500 dilution

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors