

T-47D whole cell lysate ab14899

1 References

製品の概要

製品名	T-47D whole cell lysate
特記事項	Extracts have been quality control tested by Western blot and the Electrophoretic Mobility Shift Assay (EMSA).

製品の特性

Mycoplasma free	Yes
製品の状態	Liquid
保存方法	Shipped on dry ice. Upon delivery aliquot and store at -80°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.90 Constituents: 0.75% Potassium chloride, 0.009% PMSF, 0.01% Magnesium chloride, 0.008% DTT, 0.48% HEPES, 20% Glycerol
背景	T47D cells carry receptors for a variety of steroids and calcitonin. They express mutant tumor suppressor protein p53 protein. Under normal culturing conditions, they express progesterone receptor constitutively and are responsive to estrogen. They are able to lose the estrogen receptor (ER) during long-term estrogen deprivation in vitro. Cells are therefore used as a model for studies of drug resistance to tamoxifen in patients with mutant p53 breast tumors.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors