

# Anti- Metapneumovirus Phosphoprotein antibody [hMPV 123] ab94803

## 3 References

### 製品の概要

|          |  |
|----------|--|
| 製品名      | Anti- Metapneumovirus Phosphoprotein antibody [hMPV 123]   |
| 製品の詳細    | Mouse monoclonal [hMPV 123] to Metapneumovirus Phosphoprotein  |
| 由来種      | Mouse  |
| 特異性      | Identification and location of expression of the phosphoprotein of human metapneumovirus (hMPV) of both types A and B.   |
| アプリケーション | <b>適用あり:</b> IHC-Fr, ELISA, WB, ICC/IF   |
| 種交差性     | <b>交差種:</b> Metapneumovirus  |
| 免疫原      | Tissue, cells or virus corresponding to Metapneumovirus Phosphoprotein. Human Metapneumovirus strain NCL-145   |
| 特記事項     | <p>Fusion partner: PS-NS/1-Ag4</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p> |

### 製品の特性

|        |   |
|--------|---|
| 製品の状態  | Liquid  |
| 保存方法   | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー  | Preservative: 0.02% Sodium azide<br>Constituent: 99.98% PBS   |
| 精製度    | Protein A purified  |
| ポリ/モノ  | モノクローナル   |
| クローン名  | hMPV 123  |
| アイソタイプ | IgG2a   |

## アプリケーション

**The Abpromise guarantee** Abpromise保証は、次のテスト済みアプリケーションにおけるab94803の使用に適用されます  
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

| アプリケーション | Abreviews | 特記事項   |
|----------|-----------|--|
| IHC-Fr   |           | Use at an assay dependent concentration.                                     |
| ELISA    |           | Use at an assay dependent concentration.                                     |
| WB       |           | Use at an assay dependent concentration. Predicted molecular weight: 33 kDa. |
| ICC/IF   |           | Use at an assay dependent concentration.                                     |

## ターゲット情報

**関連性** Human metapneumovirus (hMPV) is a member of the Paramyxoviridae family belonging to the Metapneumovirus genus.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

## Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

## Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors