

Anti-14-3-3 Binding Domain (phospho S) antibody ab14127

1 References

製品の概要

製品名	Anti-14-3-3 Binding Domain (phospho S) antibody
製品の詳細	Rabbit polyclonal to 14-3-3 Binding Domain (phospho S)
由来種	Rabbit
特異性	Ab14127 binds peptides and proteins containing a motif composed of phospho-Ser with proline at the +2 position and arginine or lysine at the -3 position. Antibody binding is phospho-specific and largely independent of other surrounding amino acids. The antibody weakly crossreacts with sequences containing phospho-Thr instead of phospho-Ser in this motif and with sequences containing phospho-Ser surrounded by Phe at the +1 position and Arg/Lys at the -3 position. No crossreactivity is observed with corresponding nonphosphorylated sequences or with other phospho-Thr/Ser/Tyr-containing motifs. By ELISA ab14127 recognizes a wide range of peptides containing the 14-3-3 binding motif and by 2D gel Western blot analysis it recognizes a large number of presumptive 14-3-3 binding proteins.
アプリケーション	適用あり: WB, IP, IHC-P, ELISA
種交差性	交差種: Species independent
免疫原	Synthetic peptide corresponding to 14-3-3 Binding Domain (phospho S).
特記事項	

This antibody provides a powerful tool for the discovery and characterization of potential 14-3-3 binding proteins containing this motif and for high throughput drug discovery.

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.50

	Constituents: 0.238% HEPES, 50% Glycerol, 0.87% Sodium chloride, BSA
精製度	Protein A purified
一次抗体 備考	This antibody provides a powerful tool for the discovery and characterization of potential 14-3-3 binding proteins containing this motif and for high throughput drug discovery.
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab14127の使用に適用されます**
 アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/1000. Predicted molecular weight: 30 kDa.
IP		1/50. On denatured proteins.
IHC-P		1/100.
ELISA		1/1000.

ターゲット情報

関連性 The 14-3-3 proteins are a highly conserved family of proteins involved in the regulation of cell survival, apoptosis, proliferation and checkpoint control. Biological regulation by 14-3-3 is mediated through phosphorylation-dependent protein-protein interactions. Two different phospho-Ser-containing motifs are found within nearly all known 14-3-3 binding proteins. Motif 1 (Arg/Lys and Ser at positions -3 and -2, phospho-Ser at position 0, and Pro at position +2) is found in critical regulatory proteins including Bad, cdc25c, FKHRL1, PKC and c-Raf.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors