

Anti-AC6 antibody ab14781

★★★★★ [1 Abreviews](#) [6 References](#) [画像数 1](#)

製品の概要

製品名	Anti-AC6 antibody
製品の詳細	Rabbit polyclonal to AC6
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide corresponding to Rat AC6 aa 1-100 (N terminal) conjugated to keyhole limpet haemocyanin. Database link: Q03343

 [Run BLAST with](#)

 [Run BLAST with](#)

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C or -80°C. Avoid freeze / thaw cycle.
バッファー	Preservative: 0.02% Sodium azide
精製度	Immunogen affinity purified
特記事項(精製)	The affinity-purified antibodies are purified on antigen-sepharose affinity column and supplied in antibody stabilization buffer containing preservatives with low viscosity and cryogenic properties.
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab14781の使用に適用されます


アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB	★★★★☆ (1)	1/500. Predicted molecular weight: 130 kDa.

ターゲット情報

機能	Membrane-bound, calcium-inhibitable adenylyl cyclase.
配列類似性	Belongs to the adenylyl cyclase class-4/guanylyl cyclase family. Contains 2 guanylate cyclase domains.
細胞内局在	Membrane.

画像



Anti-AC6 antibody (ab14781) at 1/500 dilution + ADCY6 Protein

Predicted band size: 130 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors