

Anti-ACADVL/VLCAD antibody [6A9AF2] ab110285

1 References [画像数 3](#)

製品の概要

製品名	Anti-ACADVL/VLCAD antibody [6A9AF2]
製品の詳細	Mouse monoclonal [6A9AF2] to ACADVL/VLCAD
由来種	Mouse
アプリケーション	適用あり: ICC/IF, IP, IHC-P
種交差性	交差種: Human
免疫原	Tissue, cells or virus. This information is considered to be commercially sensitive.
ポジティブ・コントロール	Human HDFn cells; Human liver mitochondria; Human cerebellum tissue.
特記事項	<p>This antibody clone is manufactured by Abcam. If you require a custom buffer formulation or conjugation for your experiments, please contact orders@abcam.com.</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> <p>Product was previously marketed under the MitoSciences sub-brand.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Do Not Freeze.
バッファー	<p>pH: 7.5</p> <p>Preservative: 0.02% Sodium azide</p> <p>Constituent: HEPES buffered saline</p>
精製度	Proprietary Purification
特記事項 (精製)	The antibody was produced in vitro using hybridoma grown in serum-free media and then purified by biochemical fractionation. Purity >95% by SDS-PAGE.
ポリ/モノ	モノクローナル

クローン名	6A9AF2
アイソタイプ	IgG1
軽鎖の種類	kappa

アプリケーション


The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab110285の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF		Use a concentration of 5 µg/ml.
IP		Use at an assay dependent concentration.
IHC-P		1/250. Perform heat mediated antigen retrieval - 1 min pressure cook in 1mmol EDTA pH8.

ターゲット情報


機能	Active toward esters of long-chain and very long chain fatty acids such as palmitoyl-CoA, mysritoyl-CoA and stearoyl-CoA. Can accomodate substrate acyl chain lengths as long as 24 carbons, but shows little activity for substrates of less than 12 carbons.
パスウェイ	Lipid metabolism; mitochondrial fatty acid beta-oxidation.
関連疾患	Defects in ACADVL are the cause of acyl-CoA dehydrogenase very long chain deficiency (ACADVLD) [MIM:201475]. ACADVLD is an autosomal recessive disease which leads to impaired long-chain fatty acid beta-oxidation. It is clinically heterogeneous, with three major phenotypes: a severe childhood form, with early onset, high mortality, and high incidence of cardiomyopathy; a milder childhood form, with later onset, usually with hypoketotic hypoglycemia as the main presenting feature, low mortality, and rare cardiomyopathy; and an adult form, with isolated skeletal muscle involvement, rhabdomyolysis, and myoglobinuria, usually triggered by exercise or fasting.
配列類似性	Belongs to the acyl-CoA dehydrogenase family.
細胞内局在	Mitochondrion inner membrane.

画像




Immunocytochemistry/ Immunofluorescence - Anti-ACADVL/VLCAD antibody [6A9AF2] (ab110285)

Immunocytochemistry image of ACADVL/VLCAD antibody [6A9AF2] (ab110285) stained human HDFn cells. The cells were paraformaldehyde fixed (4%, 20 min) and Triton X-100 permeabilized (0.1%, 15min). The cells were incubated with ab110285 at 5 µg/ml for 2h at room temperature or over night at 4°C. The secondary antibody was (red) Alexa Fluor® 594 goat anti-mouse IgG (H+L) used at a 1/1000 dilution for 1h. 10% Goat serum was used as the blocking agent for all blocking steps. DAPI was used to stain the cell nuclei (blue). Target protein locates mainly in mitochondria.



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-ACADVL/VLCAD antibody [6A9AF2] (ab110285)

ab110285, at 1/250 dilution, staining ACADVL/VLCAD in formalin-fixed, paraffin-embedded human cerebellum tissue. Immunoactivity is most intense in neuronal cell bodies, most notably in the large Purkinje cells.



ACADVL/VLCAD immunocaptured from 0.75mg of human liver mitochondria lysate /10ul antibody conjugated beads using ab110285.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors