

Anti-Adenine Nucleotide Translocator 2/ANT 2 antibody ab192410

画像数 1

製品の概要

製品名	Anti-Adenine Nucleotide Translocator 2/ANT 2 antibody
製品の詳細	Rabbit polyclonal to Adenine Nucleotide Translocator 2/ANT 2
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Rat, Human
免疫原	Synthetic peptide within Human Adenine Nucleotide Translocator 2/ANT 2. The exact sequence is proprietary. Database link: P05141
ポジティブ・コントロール	A549 and PC12 whole cell lysates.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.05% Sodium azide Constituent: 99% PBS
精製度	Immunogen affinity purified
特記事項 (精製)	ab192410 was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen and the purity is > 95% (by SDS-PAGE).
ポリ/モノ	ポリクローナル

アプリケーション


The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab192410の使用に適用されず**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/500. Predicted molecular weight: 33 kDa.

ターゲット情報

機能	Catalyzes the exchange of ADP and ATP across the mitochondrial inner membrane. As part of the mitotic spindle-associated MMXD complex it may play a role in chromosome segregation.
配列類似性	Belongs to the mitochondrial carrier family. Contains 3 Solcar repeats.
細胞内局在	Mitochondrion inner membrane.

画像



All lanes : Anti-Adenine Nucleotide Translocator 2/ANT 2 antibody (ab192410)

Lane 1 : A549 whole cell lysate

Lane 2 : PC12 whole cell lysate

Predicted band size: 33 kDa

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors