

Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker ab52866

リコンビナント **RabMAb**

★★★★★ **29 Abreviews** **208 References** 画像数 18

製品の概要

製品名	Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker
製品の詳細	Rabbit monoclonal [EP1332Y] to alpha Tubulin - Microtubule Marker
由来種	Rabbit
特異性	This antibody is expected to recognise most alpha tubulin proteins and not only TUBA4A.
アプリケーション	適用あり: ICC/IF, Flow Cyt (Intra), WB, IHC-P
種交差性	交差種: Mouse, Rat, Human, Pig, Drosophila melanogaster
免疫原	Synthetic peptide within Human alpha Tubulin aa 1-100 (N terminal). The exact sequence is proprietary. Database link: P68366
ポジティブ・コントロール	WB- HeLa, HEK-293, HepG2, Caco2, NIH/3T3, PC-12, RAW 264.7, PC-12, C6 Jurkat and HEK-293T whole cell lysates; human fetal kidney lysate; Mouse and rat brain lysate; Pig skeletal muscle lysates; IHC-P: Pig kidney tissue; rat kidney tissue; mouse kidney tissue; human breast cancer and stomach tissue; IHC-Fr: Rat kidney tubule tissue; Flow Cyt (intra): HepG2 cells; ICC/IF: HUVEC, HeLa and 293 cells.
特記事項	This product is a recombinant monoclonal antibody, which offers several advantages including: <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production For more information see here . Our RabMAb [®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents .

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot. Store at -20°C. Avoid freeze / thaw cycle.
バッファー	pH: 7.20 Preservative: 0.01% Sodium azide

	Constituents: 59% PBS, 40% Glycerol (glycerin, glycerine), 0.5% BSA
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EP1332Y
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab52866の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
ICC/IF	★★★★★ (7)	1/250 - 1/500.
Flow Cyt (Intra)		1/20 - 1/50. ab172730 - Rabbit monoclonal IgG, is suitable for use as an isotype control with this antibody.
WB	★★★★★ (19)	1/1000 - 1/50000. Detects a band of approximately 50 kDa (predicted molecular weight: 50 kDa).
IHC-P	★★★★★ (1)	Use at an assay dependent concentration. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

ターゲット情報

機能	Tubulin is the major constituent of microtubules. It binds two moles of GTP, one at an exchangeable site on the beta chain and one at a non-exchangeable site on the alpha chain.
配列類似性	Belongs to the tubulin family.
翻訳後修飾	<p>Some glutamate residues at the C-terminus are polyglutamylated. This modification occurs exclusively on glutamate residues and results in polyglutamate chains on the gamma-carboxyl group. Also monoglycylated but not polyglycylated due to the absence of functional TTL10 in human. Monoglycylation is mainly limited to tubulin incorporated into axonemes (cilia and flagella) whereas glutamylation is prevalent in neuronal cells, centrioles, axonemes, and the mitotic spindle. Both modifications can coexist on the same protein on adjacent residues, and lowering glycylation levels increases polyglutamylation, and reciprocally. The precise function of such modifications is still unclear but they regulate the assembly and dynamics of axonemal microtubules.</p> <p>Acetylation of alpha chains at Lys-40 stabilizes microtubules and affects affinity and processivity of microtubule motors. This modification has a role in multiple cellular functions, ranging from cell motility, cell cycle progression or cell differentiation to intracellular trafficking and signaling.</p>
細胞内局在	Cytoplasm > cytoskeleton.

画像




Immunocytochemistry/ Immunofluorescence - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)

Immunofluorescent analysis of 4% paraformaldehyde-fixed HeLa (human epithelial cell line from cervix adenocarcinoma) cells labeling alpha Tubulin with ab52866 at 1/500 dilution. The cells were permeabilised with 0.1% Triton X-100. Anti-rabbit Alexa Fluor® 488 ([ab150077](#)) at 1/400 dilution was used as the secondary antibody (green). The confocal image shows microtubules staining on HeLa cell line. The nuclear counter stain is DAPI (blue). Tubulin is detected with [ab7291](#) (anti-Tubulin mouse mAb) at 1/500 and anti-mouse AlexaFluor® 594 ([ab150120](#)) at 1/500 dilution (red).

The negative controls are as follows:

1. ab52866 at 1/500 dilution followed by anti-mouse AlexaFluor® 594 ([ab150120](#)) at 1/500 dilution.
2. [ab7291](#) (anti-Tubulin mouse mAb) at 1/500 dilution followed by anti-rabbit Alexa Fluor® 488 ([ab150077](#)) at 1/400 dilution.



Western blot - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)

All lanes : Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866) at 1/10000 dilution

Lane 1 : HeLa (Human cervix adenocarcinoma epithelial cell)

Whole cell lysates

Lane 2 : PC-12 (Rat adrenal gland pheochromocytoma) whole cell lysates

Lane 3 : NIH/3T3(Mouse embryonic fibroblast)whole cell lysates

Lysates/proteins at 15 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG H&L (HRP) ([ab97051](#)) at 1/20000 dilution

Predicted band size: 50 kDa

Observed band size: 52 kDa

Exposure time: 1 second


Blocking/Diluting buffer and concentration: 5% NFDM/TBST



Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)


Immunohistochemical analysis of paraffin-embedded Human stomach tissue labeling alpha Tubulin with ab52866 at followed by a ready to use Goat Anti-Rabbit IgG H&L (HRP). Cytoplasmic staining on human stomach. The section was incubated with **ab229902** for 30 mins at room temperature. The immunostaining was performed on a Leica Biosystems BOND® RX instrument. Counterstained with Hematoxylin. Heat mediated antigen retrieval using **ab93684** (Tris/EDTA buffer, pH 9.0).

Secondary antibody only control: Used PBS instead of primary antibody, secondary antibody is a ready to use Goat Anti-Rabbit IgG H&L (HRP).



Flow Cytometry (Intracellular) - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)

Intracellular flow cytometric analysis of 4% paraformaldehyde fixed 90% methanol permeabilized HepG2 (Human hepatocellular carcinoma epithelial cell) cells labelling alpha Tubulin with ab52866 at 1/2000 dilution (Red) compared with a Rabbit monoclonal IgG (**ab172730**) isotype control (black) and an unlabelled control (cells without incubation with primary antibody and secondary antibody) (Blue). A Goat anti rabbit IgG (Alexa Fluor® 488, **ab150077**) at 1/2000 dilution was used as the secondary antibody.



All lanes : Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866) at 1/5000 dilution

Lane 1 : Mouse brain lysates

Lane 2 : C6 (Rat glial tumor cell line) whole cell lysates

Lane 3 : RAW 264.7 (Mouse macrophage cell line transformed with Abelson murine leukemia virus) whole cell lysates

Lane 4 : PC-12 (Rat adrenal gland pheochromocytoma cell line) whole cell lysates


Lane 5 : NIH/3T3 (Mouse embryonic fibroblast cell line) whole cell lysates

Lysates/proteins at 10 µg per lane.

Secondary

All lanes : Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at 1/1000 dilution

Predicted band size: 50 kDa




Immunohistochemistry analysis of paraffin-embedded Pig kidney tissue labeling alpha Tubulin with ab52866 at a 1/1000 dilution. Cytoplasmic staining on Pig kidney tubule and weak on glomerulus shown. **Anti-Rabbit HRP (ab97051)** used at a 1/100 dilution. Counter stained with Hematoxylin.

Inset image: negative control obtained using PBS instead of ab52866, secondary antibody is **Anti-Rabbit HRP (ab97051)** at 1/100 dilution.


Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)



Flow Cytometry (Intracellular) - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)

Intracellular Flow Cytometry analysis of 2% paraformaldehyde fixed HepG2 (human liver hepatocellular carcinoma cell line) cells labeling alpha Tubulin with ab52866 at 1/130 dilution (red line). Secondary antibody used is a goat anti rabbit IgG (FITC) at 1/150 dilution. The isotype control is rabbit monoclonal IgG (black line). The unlabeled control is cells without incubation with primary and secondary antibodies (blue line).




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)

Immunohistochemistry analysis of paraffin-embedded Rat kidney tissue labeling alpha Tubulin with ab52866 at a 1/1000 dilution. Cytoplasmic staining on Rat kidney tubule and weak on glomerulus shown. Secondary antibody **Anti-Rabbit HRP (ab97051)** used at a 1/500 dilution. Counter stained with Hematoxylin.

Inset image: negative control obtained using PBS instead of ab52866, secondary antibody is **Anti-Rabbit HRP (ab97051)** at 1/500 dilution.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.



Immunocytochemistry/ Immunofluorescence - Anti-
alpha Tubulin antibody [EP1332Y] - Microtubule
Marker (ab52866)


Hotte, K. et al Biol Open. 2017 Jun 15;6(6):723-731. doi:
10.1242/bio.022335 Reproduced under the Creative
Commons license
<http://creativecommons.org/licenses/by/3.0>

AJAP1 co-localizes with microtubules in HUVECs

The association of AJAP1 with microtubules in HUVECs is lost upon microtubule destruction. Treatment with 12.5 μ M nocodazole for 24 h shows destruction of the microtubule network and loss of AJAP1 tubular localization. For a negative control, HUVECs are treated with DMSO for 24 h. Cell nuclei were counterstained with DAPI (cyan). Microscope: Zeiss LSM 780; objective lens: 63 \times /1.40 oil; scale bar: 25 μ m.

Incubated overnight at 4°C with ab52866.

(From Figure 3E of Hotte et al)



Western blot - Anti-alpha Tubulin antibody
[EP1332Y] - Microtubule Marker (ab52866)

All lanes : Anti-alpha Tubulin antibody [EP1332Y] - Microtubule
Marker (ab52866) at 1/20000 dilution

Lane 1 : HeLa (human epithelial cell line from cervix
adenocarcinoma) whole cell lysate

Lane 2 : HepG2 (human liver hepatocellular carcinoma cell line)
whole cell lysate

Lane 3 : Jurkat (human T cell leukemia cell line from peripheral
blood) whole cell lysate


Lane 4 : HEK-293T (human epithelial cell line from embryonic
kidney transformed with large T antigen) whole cell lysate

Lysates/proteins at 20 μ g per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at
1/1000 dilution

Predicted band size: 50 kDa




Western blot - Anti-alpha Tubulin antibody
[EP1332Y] - Microtubule Marker (ab52866)

Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker
(ab52866) at 1/50000 dilution + Rat brain lysates at 10 µg

Secondary

Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at
1/1000 dilution


Predicted band size: 50 kDa



Immunocytochemistry/ Immunofluorescence - Anti-
alpha Tubulin antibody [EP1332Y] - Microtubule
Marker (ab52866)

This image is courtesy of an anonymous Abreview

ab52866 staining alpha Tubulin in 293 Human embryonic kidney cells by ICC/IF (Immunocytochemistry/immunofluorescence). Cells were fixed with paraformaldehyde and blocked with 10% serum for 2 hours at 23°C. Samples were incubated with primary antibody (1/200 in 0.5% saponin) for 2 hours at 23°C. An Alexa Fluor®555-conjugated Goat anti-rabbit IgG polyclonal (1/500) was used as the secondary antibody. Nuclei were counterstained with DAPI.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)

Immunohistochemistry analysis of paraffin-embedded Mouse kidney tissue labeling alpha Tubulin with ab52866 at a 1/1000 dilution. Cytoplasmic staining on Mouse kidney tubule shown. Secondary antibody **Anti-Rabbit HRP (ab97051)** used at a 1/500 dilution. Counter stained with Hematoxylin.

Inset image: negative control obtained using PBS instead of ab52866, secondary antibody is **Anti-Rabbit HRP (ab97051)** at 1/500 dilution.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.




Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker (ab52866)

Immunohistochemistry analysis of paraffin-embedded Human breast cancer labeling alpha Tubulin with ab52866 at a 1/1000 dilution. Cytoplasmic staining on cancer cells shown. Secondary antibody **ab97051** Goat Anti-Rabbit IgG H&L (HRP) used at a 1/500 dilution. Counter stained with Hematoxylin.

Inset image: negative control obtained using PBS instead of ab52866, secondary antibody is Anti-Rabbit HRP (**ab97051**) at 1/500 dilution.

Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.




Western blot - Anti-alpha Tubulin antibody
[EP1332Y] - Microtubule Marker (ab52866)

Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker
(ab52866) at 1/50000 dilution + Human fetal kidney lysates at 10 µg

Secondary

Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at
1/1000 dilution

Predicted band size: 50 kDa



Western blot - Anti-alpha Tubulin antibody
[EP1332Y] - Microtubule Marker (ab52866)

All lanes : Anti-alpha Tubulin antibody [EP1332Y] - Microtubule
Marker (ab52866) at 1/1000 dilution

Lane 1 : HeLa (Human epithelial carcinoma cell line) Whole Cell
Lysate

Lane 2 : HEK293 (Human embryonic kidney cell line) Whole Cell
Lysate

Lane 3 : HepG2 (Human hepatocellular liver carcinoma cell line)
Whole Cell Lysate

Lane 4 : Caco 2 (Human colonic carcinoma cell line) Whole Cell
Lysate

Lane 5 : NIH/3T3 (Mouse embryonic fibroblast cell line) Whole Cell
Lysate

Lane 6 : PC-12 (Rat adrenal pheochromocytoma cell line) Whole
Cell Lysate

Lysates/proteins at 20 µg per lane.

Secondary


All lanes : Goat Anti-Rabbit IgG H&L (Alexa Fluor® 790)
(**ab175781**) at 1/10000 dilution

Predicted band size: 50 kDa

Observed band size: 52 kDa

This blot was produced using a 4-12% Bis-tris gel under the MOPS
buffer system. The gel was run at 200V for 50 minutes before being
transferred onto a Nitrocellulose membrane at 30V for 70 minutes.

The membrane was then blocked for an hour using Licor blocking buffer before being incubated with ab52866 overnight at 4°C. Antibody binding was detected using **Anti-Rabbit Alexa Fluor® 790 (ab175781)** at a 1:10,000 dilution for 1hr at room temperature and then imaged using the Licor Odyssey CLx.



Western blot - Anti-alpha Tubulin antibody
[EP1332Y] - Microtubule Marker (ab52866)

Anti-alpha Tubulin antibody [EP1332Y] - Microtubule Marker
(ab52866) at 1/5000 dilution + Pig skeletal muscle lysates at 20 µg

Secondary

Anti-Rabbit IgG (HRP), specific to the non-reduced form of IgG at
1/1000 dilution

Predicted band size: 50 kDa

Why choose a recombinant antibody?



Anti-alpha Tubulin antibody [EP1332Y] - Microtubule
Marker (ab52866)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors