


Anti-Amphiphysin antibody ab96508

面格新?

医薬用外劇物

製品の概要

製品名	Anti-Amphiphysin antibody
製品の詳細	Rabbit polyclonal to Amphiphysin
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Cow 
免疫原	Recombinant fragment, corresponding to a region within the N terminal amino acids 1-164 of Human Amphiphysin.
ポジティブ・コントロール	293T, A431, H1299, HeLa S3, HepG2, MOLT4, Raji cell lysates NCIN87 xenograft (IHC)
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 1.21% Tris, 0.75% Glycine, 10% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab96508の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 76 kDa.
IHC-P		1/100 - 1/500. Antigen retrieval is not essential but may optimise staining.

ターゲット情報

機能	May participate in mechanisms of regulated exocytosis in synapses and certain endocrine cell types. May control the properties of the membrane associated cytoskeleton.
組織特異性	Neurons, certain endocrine cell types and spermatocytes.
配列類似性	Contains 1 BAR domain. Contains 1 SH3 domain.
細胞内局在	Cytoplasmic vesicle > secretory vesicle > synaptic vesicle membrane. Cytoplasm > cytoskeleton.


画像



Anti-Amphiphysin antibody (ab96508) at 1/1000 dilution + 293T whole cell lysate at 30 µg

Predicted band size: 76 kDa

7.5% SDS PAGE



Immunohistochemical analysis of paraffin-embedded NCIN87 xenograft, using ab96508 at 1:100 dilution.

Immunohistochemistry (Formalin/PFA-fixed paraffin-embedded sections) - Anti-Amphiphysin antibody (ab96508)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors