

Anti-Angiopoietin 1 antibody [EPR2888(N)] ab183701

リコンビナント **RabMAb**

16 References 画像数 3

製品の概要

製品名	Anti-Angiopoietin 1 antibody [EPR2888(N)]
製品の詳細	Rabbit monoclonal [EPR2888(N)] to Angiopoietin 1
由来種	Rabbit
アプリケーション	適用あり: WB
種交差性	交差種: Mouse, Rat, Human
免疫原	Synthetic peptide. This information is proprietary to Abcam and/or its suppliers.
ポジティブ・コントロール	WB: U87-MG, HeLa, C6, RAW 264.7, K562, and ECV-304 cell lysates.
特記事項	<p>This product is a recombinant monoclonal antibody, which offers several advantages including:</p> <ul style="list-style-type: none"> - High batch-to-batch consistency and reproducibility - Improved sensitivity and specificity - Long-term security of supply - Animal-free production <p>For more information see here.</p> <p>Our RabMAb[®] technology is a patented hybridoma-based technology for making rabbit monoclonal antibodies. For details on our patents, please refer to RabMAb[®] patents.</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C short term (1-2 weeks). Upon delivery aliquot. Store at -20°C long term. Avoid freeze / thaw cycle.
バッファー	<p>pH: 7.2</p> <p>Preservative: 0.01% Sodium azide</p> <p>Constituents: 40% Glycerol, 0.05% BSA, 59% PBS</p>
精製度	Protein A purified
ポリ/モノ	モノクローナル
クローン名	EPR2888(N)
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab183701の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
WB		1/10000 - 1/50000. Detects a band of approximately 70 kDa (predicted molecular weight: 57 kDa).

ターゲット情報

機能

Binds and activates TIE2 receptor by inducing its tyrosine phosphorylation. Implicated in endothelial developmental processes later and distinct from that of VEGF. Appears to play a crucial role in mediating reciprocal interactions between the endothelium and surrounding matrix and mesenchyme. Mediates blood vessel maturation/stability. It may play an important role in the heart early development.

配列類似性

Contains 1 fibrinogen C-terminal domain.


翻訳後修飾

Glycosylated.

細胞内局在

Secreted.

画像



All lanes : Anti-Angiopoietin 1 antibody [EPR288(N)] (ab183701) at 1/10000 dilution (Purified)

Lane 1 : U-87 MG (Human glioblastoma-astrocytoma epithelial cell) whole cell lysates

Lane 2 : C6 (Rat glial tumor glial cell) whole cell lysates

Lane 3 : RAW 264.7 (Mouse Abelson murine leukemia virus-induced tumor macrophage) whole cell lysates

Lane 4 : Human brain lysates

Lane 5 : Mouse brain lysates

Lane 6 : Rat brain lysates

Lane 7 : Human heart lysates

Lane 8 : Mouse heart lysates

Lysates/proteins at 20 µg per lane.


Secondary

All lanes : Goat Anti-Rabbit IgG H&L (HRP) (**ab97051**) at 1/20000 dilution

Predicted band size: 57 kDa

Observed band size: 70 kDa

This antibody is unsuitable to detect tissue lysate samples.



All lanes : Anti-Angiopoietin 1 antibody [EPR2888(N)] (ab183701)
at 1/50000 dilution

Lane 1 : U87-MG cell lysate

Lane 2 : HeLa cell lysate

Lane 3 : K562 cell lysate

Lane 4 : ECV-304 cell lysate

Lysates/proteins at 20 µg per lane.

Secondary

All lanes : Goat Anti-Rabbit IgG, (H+L), Peroxidase conjugated at
1/1000 dilution

Predicted band size: 57 kDa

Why choose a recombinant antibody?

Research with confidence
Consistent and reproducible results

Long-term and scalable supply
Recombinant technology

Success from the first experiment
Confirmed specificity

Ethical standards compliant
Animal-free production

Anti-Angiopoietin 1 antibody [EPR2888(N)]
(ab183701)

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery

- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors