


Anti-Annexin-8/ANXA8 antibody ab111693

1 References [画像数 3](#)

医薬用外劇物

製品の概要

製品名	Anti-Annexin-8/ANXA8 antibody
製品の詳細	Rabbit polyclonal to Annexin-8/ANXA8
由来種	Rabbit
アプリケーション	適用あり: WB, IHC-P, ICC/IF
種交差性	交差種: Human 交差が予測される動物種: Mouse, Rat, Rabbit 
免疫原	Recombinant fragment corresponding to Human Annexin-8/ANXA8 aa 1-188. (AAH04376)
ポジティブ・コントロール	A549 and HeLa whole cell lysates; HeLa cells; U87 xenograft tissue.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles.
バッファー	pH: 7.00 Preservative: 0.01% Thimerosal (merthiolate) Constituents: 78.99% PBS, 1% BSA, 20% Glycerol (glycerin, glycerine)
精製度	Immunogen affinity purified
ポリ/モノ	ポリクローナル
アイソタイプ	IgG

アプリケーション


The Abpromise guarantee **Abpromise保証は、**次のテスト済みアプリケーションにおけるab111693の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
WB		1/500 - 1/3000. Predicted molecular weight: 37 kDa.
IHC-P		1/100 - 1/500. Perform heat mediated antigen retrieval with citrate buffer pH 6 before commencing with IHC staining protocol.
ICC/IF		1/100 - 1/500.

ターゲット情報


機能	This protein is an anticoagulant protein that acts as an indirect inhibitor of the thromboplastin-specific complex, which is involved in the blood coagulation cascade.
配列類似性	Belongs to the annexin family. Contains 4 annexin repeats.
ドメイン	A pair of annexin repeats may form one binding site for calcium and phospholipid.

画像



ab111693 at 1/200 dilution staining Annexin-8/ANXA8 in methanol-fixed HeLa cells by Immunofluorescence. Lower image shows cells co-stained with Hoechst 33342.

Immunocytochemistry/ Immunofluorescence - Anti-Annexin-8/ANXA8 antibody (ab111693)




Western blot - Anti-Annexin-8/ANXA8 antibody
(ab111693)

Anti-Annexin-8/ANXA8 antibody (ab111693) at 1/1000 dilution +
A549 whole cell lysate at 30 µg

Predicted band size: 37 kDa

10% SDS-PAGE.



Immunohistochemistry (Formalin/PFA-fixed paraffin-
embedded sections) - Anti-Annexin-8/ANXA8
antibody (ab111693)

ab111693 at 1/500 dilution staining Annexin-8/ANXA8 in Formalin-
fixed, paraffin-embedded U87 xenograft tissue by
Immunohistochemistry.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise,

please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors