

APC Anti-CD90 / Thy1 antibody [G7] ab25322

2 References [画像数 1](#)

製品の概要

製品名	APC Anti-CD90 / Thy1 antibody [G7]
製品の詳細	APC Rat monoclonal [G7] to CD90 / Thy1
由来種	Rat
標識	APC. Ex: 645nm, Em: 660nm
特異性	ab25322 recognises Thy1.1 and Thy1.2 alloantigens.
アプリケーション	適用あり: Flow Cyt
種交差性	交差種: Mouse
免疫原	Tissue, cells or virus corresponding to Mouse CD90/ Thy1. Mouse T cell hybridoma C6/G8
エピトープ	ab25322 reacts with an epitope located in the region A of CD90/Thy1.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	<p>pH: 7.3</p> <p>Preservative: 0.09% Sodium azide</p> <p>Constituents: PBS, 16% Sucrose</p> <p>Stabilising agent.</p>
精製度	Affinity purified
ポリ/モノ	モノクローナル
クローン名	G7
アイソタイプ	IgG2c

アプリケーション


The Abpromise guarantee Abpromise保証は、次のテスト済みアプリケーションにおけるab25322の使用に適用されます
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		Use 0.2µg for 10 ⁶ cells.

ターゲット情報

機能	May play a role in cell-cell or cell-ligand interactions during synaptogenesis and other events in the brain.
配列類似性	Contains 1 Ig-like V-type (immunoglobulin-like) domain.
細胞内局在	Cell membrane.

画像



Flow cytometry analysis of BALB/c mouse splenocytes stained with Rat Anti-Mouse CD90-APC ab25322 (0.01 µg/10⁶cells) and Rat Anti-Mouse CD3 epsilon-FITC [ab24948](#).

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors