

# APC Mouse IgG2b [PLPV219] - Isotype Control ab91534

## 製品の概要

製品名	APC Mouse IgG2b [PLPV219] - Isotype Control
標識	APC. Ex: 645nm, Em: 660nm
特異性	ab91534 was produced against a synthetic hapten, which is normally not present in Humans or animals. Therefore, ab91534 does not react with any antigen of Human or animal cells.
アプリケーション	<b>適用あり:</b> Flow Cyt
免疫原	Vegetal protein
特記事項	<p>25µg is sufficient for 50 tests</p> <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&amp;As</p>

## 製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Store at +4°C. Store In the Dark.
バッファー	<p>pH: 7.20</p> <p>Preservative: 0.1% Sodium azide</p> <p>Constituents: PBS, 1% BSA</p>
精製度	Protein G purified
ポリ/モノ	モノクローナル
クローン名	PLPV219
アイソタイプ	IgG2b

## アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab91534の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

アプリケーション	Abreviews	特記事項
Flow Cyt		

#### 追加情報

Flow Cyt: Use 20µl for 10<sup>6</sup> cells.  
ab91534 is a negative isotype control

Not yet tested in other applications.  
Optimal dilutions/concentrations should be determined by the end user.

**Please note:** All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

#### Our Abpromise to you: Quality guaranteed and expert technical support

---

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

#### Terms and conditions

---

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors