


Anti-APG8A/ATG8A antibody ab77003

★★★★★ **1 Abreviews** **30 References**

製品の概要

| | |
|-------|---|
| 製品名 | Anti-APG8A/ATG8A antibody |
| 製品の詳細 | Rabbit polyclonal to APG8A/ATG8A |
| 由来種 | Rabbit |
| 種交差性 | 交差種: Arabidopsis thaliana, Corn 交差が予測される動物種: Human  |
| 免疫原 | Recombinant full length protein corresponding to APG8A/ATG8A. |
| 特記事項 | <p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p> |

製品の特性

| | |
|----------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C or -80°C. Avoid repeated freeze / thaw cycles. |
| バッファー | Constituents: 2.42% Tris, 1.5% Glycine |
| 精製度 | IgG fraction |
| 特記事項(精製) | This antibody was purified over AFFI-GEL 10 beads binds with rATG8 protein. |
| ポリ/モノ | ポリクローナル |
| アイソタイプ | IgG |

ターゲット情報

| | |
|-----|--|
| 関連性 | APG8A is involved in cytoplasm to vacuole transport (Cvt) vesicles and autophagosomes formation. It may mediate the delivery of the vesicles and autophagosomes to the vacuole via the microtubule cytoskeleton. |
|-----|--|

細胞内局在

Cytoplasmic vesicle › cvt vesicle membrane; Lipid-anchor. Cytoplasmic vesicle › autophagosome membrane; Lipid-anchor. Vacuole membrane; Lipid-anchor. Note= Membrane-associated through a lipid anchor.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors