

Anti-Bacillus anthracis protective antigen antibody [PA6] ab69485

製品の概要

| | |
|----------|--|
| 製品名 | Anti-Bacillus anthracis protective antigen antibody [PA6] |
| 製品の詳細 | Mouse monoclonal [PA6] to Bacillus anthracis protective antigen |
| 由来種 | Mouse |
| アプリケーション | 適用あり: ELISA, WB, IP |
| 種交差性 | 交差種: Bacillus anthracis |
| 免疫原 | Full length protein corresponding to Bacillus anthracis protective antigen. Database link: P13423 |

特記事項

The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.

If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As

製品の特性

| | |
|--------|---|
| 製品の状態 | Liquid |
| 保存方法 | Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles. |
| バッファー | Preservative: 0.01% Sodium azide Constituent: PBS |
| 精製度 | Protein G purified |
| ポリ/モノ | モノクローナル |
| クローン名 | PA6 |
| アイソタイプ | IgG1 |
| 軽鎖の種類 | kappa |

アプリケーション

The Abpromise guarantee

Abpromise保証は、次のテスト済みアプリケーションにおけるab69485の使用に適用されます

アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご確認ください。

| アプリケーション | Abreviews | 特記事項 |
|----------|-----------|---|
| ELISA | | Use at an assay dependent dilution. |
| WB | | Use at an assay dependent dilution. Predicted molecular weight: 86 kDa. |
| IP | | Use at an assay dependent dilution. |

ターゲット情報

関連性 Anthrax toxins are composed of three distinct proteins, a protective antigen (PA), a lethal factor (LF) and an edema factor (EF). None of these proteins are toxic by themselves and several studies indicate that the anthrax toxin has the familiar A-B enzymatic binding structure, with PA acting as the binding domain and EF and/or LF acting as the active fragments.

細胞内局在 Secreted

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours

- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors