

Anti-Bacillus thuringiensis CRY1Ab Toxin antibody ab51586

4 References

製品の概要

製品名	Anti-Bacillus thuringiensis CRY1Ab Toxin antibody
製品の詳細	Rabbit polyclonal to Bacillus thuringiensis CRY1Ab Toxin
由来種	Rabbit
アプリケーション	適用あり: IP, WB, ELISA
種交差性	交差種: Bacillus thuringiensis
免疫原	Full length native protein (purified) corresponding to Bacillus thuringiensis CRY1Ab Toxin.
特記事項	<p>The Life Science industry has been in the grips of a reproducibility crisis for a number of years. Abcam is leading the way in addressing this with our range of recombinant monoclonal antibodies and knockout edited cell lines for gold-standard validation. Please check that this product meets your needs before purchasing.</p> <p>If you have any questions, special requirements or concerns, please send us an inquiry and/or contact our Support team ahead of purchase. Recommended alternatives for this product can be found below, along with publications, customer reviews and Q&As</p>

製品の特性

製品の状態	Liquid
保存方法	Shipped at 4°C. Upon delivery aliquot and store at -20°C. Avoid freeze / thaw cycles.
バッファー	pH: 7.30 Constituent: PBS
精製度	Protein A purified
ポリモノ	ポリクローナル
アイソタイプ	IgG

アプリケーション

The Abpromise guarantee **Abpromise保証は、次のテスト済みアプリケーションにおけるab51586の使用に適用されます**
アプリケーションノートには、推奨の開始希釈率がありますが、適切な希釈率につきましてはご検討ください。

アプリケーション	Abreviews	特記事項
IP		Use at an assay dependent concentration.
WB		Use at an assay dependent concentration. Predicted molecular weight: 131 kDa.
ELISA		1/500 - 1/800.

ターゲット情報

関連性 CRY1A belongs to the delta endotoxin family and is one of the most studied insecticidal proteins produced by *Bacillus thuringiensis*. The crystal protein is produced during sporulation and is accumulated both as an inclusion and as part of the spore coat. The toxic segment of the protein is located in the N-terminus. Promotes colloid osmotic lysis by binding to the midgut epithelial cells of many lepidopteran larvae.

Please note: All products are "FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES"

Our Abpromise to you: Quality guaranteed and expert technical support

- Replacement or refund for products not performing as stated on the datasheet
- Valid for 12 months from date of delivery
- Response to your inquiry within 24 hours
- We provide support in Chinese, English, French, German, Japanese and Spanish
- Extensive multi-media technical resources to help you
- We investigate all quality concerns to ensure our products perform to the highest standards

If the product does not perform as described on this datasheet, we will offer a refund or replacement. For full details of the Abpromise, please visit <https://www.abcam.co.jp/abpromise> or contact our technical team.

Terms and conditions

- Guarantee only valid for products bought direct from Abcam or one of our authorized distributors